

Contact information:

3Shape Implant Studio



Drilling Protocol


Order Details

Patient Name: Sagled Berriche

Client Order Reference: 409450643_20210611_1103_22_067

Creation Date: 11/06/2021 11:51:27

Created by: Lallouette La Perle, labolaperle@gmail.com
(b2b2b2b2-b2b2-b2b2-8c35-a843012c8c4b)



*1>Note: Please check the required final drill length against the suggested minimum drill length for your chosen implants. The final drill length depends on the chosen sleeve type and instructions for use for the relevant surgical protocol defined by your implant choice. E.g., for fully guided sleeves additional values might need to be added like a potentially used spoon or extended drill depth.

Implant information	
Implant position (FDI)	25
Manufacturer	Biotech Dental
Type	Biotech Dental KONTACT Internal 4.2x12.0
Order number	K4212
Length, mm	12
Diameter (Ø), mm	4.2
Color	Blue
Sleeve information	
Name	Biotech Dental KMTW
Type	Fully guided
Order number	KMTW
Offset, mm	9
Drill information	
Minimum drill length	21

Limitation of Liability:

This instruction incorporates a custom document that is based on a surgical plan proposed by the surgeon before operation. The surgeon, therefore, takes full medical responsibility for the design and the application of the surgical guide, the intended used surgical tray kit, implants and sleeves – all as specified on the order form received by the supplier. The custom document shall be considered as an addition to all other documents sent with and pertaining to the case, and it does not replace any of those other documents.